


# Common Fixation Pins

Stainless steel material meets ASTM F138 and ISO 5832-1 standards. Available partially threaded or smooth. Also available pointed one-end only with plain or radius end.

Customized to your specific requirements.



Shaft Diameter (A)		Overall Length (B)		Thread Size (C)
inches	millimeters	inches	millimeters	
0.063	1.6	7.0	177.80	0-80
0.079	2.0	7.0	177.80	2-56
0.094	2.4	8.0	203.20	3-48
0.110	2.8	8.0	203.20	4-40
0.125	3.2	9.0	228.60	5-40
0.141	3.6	9.0	228.60	6-32
0.157	4.0	10.0	254.00	8-32
0.188	4.8	12.0	304.80	10-24
0.250	6.4	12.0	304.80	1/4-20