




Wound Drainage Trocars

Common wound drainage trocars are manufactured continuously and available from stock.



Part Number	Connection Type (A)	Shaft Diameter (B)	Overall Length (C)	Connection Length (D)	Connection Diameter (E)
102-001B-0236-152-00	Barbed	3/32 in. 2.381 mm.	6.0 in. 152.400 mm.	0.266 in. 6.756 mm.	.070 in. 1.778 mm.
102-001B-0318-152-00	Barbed	1/8 in. 3.175 mm.	6.0 in. 152.400 mm.	0.266 in. 6.756 mm.	0.092 in. 2.337 mm.
102-001B-0397-152-00	Barbed	5/32 in. 3.97 mm.	6.0 in. 152.400 mm.	0.266 in. 6.756 mm.	0.110 in. 2.794 mm.
102-001B-0476-152-00	Barbed	3/16 in. 4.763 mm.	6.0 in. 152.400 mm.	0.440 in. 11.176 mm.	0.159 in. 4.039 mm.
102-001B-0476-152-10	Small Barbed	3/16 in. 4.763 mm.	6.0 in. 152.400 mm.	0.427 in. 10.846 mm.	0.123 in. 3.124 mm.
102-001B-0635-152-00	Barbed	1/4 in. 6.350 mm.	6.0 in. 152.400 mm.	0.440 in. 11.176 mm.	0.159 in. 4.039 mm.
102-001B-0635-152-10	Special Barbed	1/4 in. 6.350 mm.	6.0 in. 152.400 mm.	0.430 in. 10.92 mm.	0.190 in. 4.83 mm.
002-0018-0635-152-20	Large Barbed	1/4 in. 6.350 mm.	6.0 in. 152.400 mm.	0.571 in. 14.50 mm.	0.225 in. 5.715 mm.
103-001C-0318-152-00	Threaded	1/8 in. 3.175 mm.	6.0 in. 152.400 mm.	0.343 in. 8.712 mm.	0.087 in. 2.210 mm.
103-001C-0476-152-00	Threaded	3/16 in. 4.763 mm.	6.0 in. 152.400 mm.	0.562 in. 14.275 mm.	0.147 in. 3.734 mm.
103-001C-0635-152-00	Threaded	1/4 in. 6.350 mm.	6.0 in. 152.400 mm.	0.562 in. 14.275 mm.	0.147 in. 3.734 mm.