



Common Handle Types Manufactured



Handles



Handle Formed



Sphere/Ball
Plastic/Metal


Common Distal Point Types Manufactured




Truncated Cone



Drilled Cone




Flat



Bullet




Trocar Point




Screw Driver




Chisel with Back Grind




Reducer



Pencil




3 Facet Point



4 Facet Point


Common Proximal End Types Manufactured




Plain




Tapered



Drilled




Barb



Reduced Diameter




Diamond Knurl




Taper




Threaded



Reverse Taper




Bent




Double Reduction




Straight Knurl




Bullet



Blunt



Double Angle



Double Flat